

EndoSafe

Safer for the patient


- Designed for maximum patient safety
- No need for closed insufflation; lower or no risk of CO₂ embolism
- No risk of saphenous vein damage during vessel harvesting of the collateral vessels due to the special design of the manipulator.



EndoEasy

Easier to use

- The Endoperfect system is very easy and simple to use
- The Spacer device is designed to house the endoscopic instruments inside the operative channel which acts as a guide without the need for CO₂ or manually lifting the instruments
- Optic housing channel separated from the spacer
- Smoke suction channel and lens washing channel to ensure the best possible visibility of the surgical field
- 30° optics: decreased risk of soiling when compared to a 0° optic



EndoLogic

Higher efficiency at a lower cost

- Reduced costs:
 - no need for specialized instrumentation
 - no need for expensive mono use optics or bipolar scalpels
- Versatile: compatible with almost all currently available endoscopic devices